

# SOP-CAL-018 Revision Recommendation — DO Probe Calibration

AMGEN — DEVIATION TREND INTELLIGENCE

SOP Revision Recommendation

Generated by: Pace AI Agent

Date: 2026-02-18

Priority: HIGH

TARGET SOP: SOP-CAL-018 Rev 4 (DO Probe Calibration Procedure)

Applicable Sites: ATO, IRE, PR

## FINDING

7 dissolved oxygen excursion NCs across 3 sites in 4 weeks, all occurring on Days 12-16 of culture (peak cell density phase). All batches showed titer >10% above historical average (4.8 g/L vs 4.2 g/L). Current SOP-CAL-018 calibration was validated for titer up to 4.5 g/L.

## ROOT CAUSE

Media optimization has increased average titer from 4.2 to 4.8 g/L. Higher titer = higher cell density = higher oxygen demand. DO probes calibrated per current SOP are accurate but response is insufficient for the increased demand profile. The probes are operating beyond their validated range.

## HISTORICAL PRECEDENT

In Q3 2024, an identical pattern (5 minor DO NCs over 6 weeks) preceded 2 batch failures at ATO, costing \$3.2M. Current trajectory (7 NCs in 4 weeks) is more aggressive.

## RECOMMENDED REVISIONS

1. IMMEDIATE: Increase DO probe calibration frequency from quarterly to monthly
2. SHORT-TERM: Expand validation range in SOP-CAL-018 Rev 5 to cover titer up to 5.5 g/L
3. LONG-TERM: Evaluate optical DO sensors (e.g., Hamilton VisiFerm) for high-titer applications

Change control to be initiated in TrackWise.